

PPT2 Antibody, Biotin conjugated



PACO64486

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Removes thioester-linked fatty acyl groups from various substrates including S-palmitoyl-CoA. Has the highest S-thioesterase activity for the acyl groups palmitic and myristic acid, followed by other short- and long-chain acyl substrates. However, because of structural constraints, is unable to remove palmitate from peptides or proteins.

Gene ID:

PPT2

Uniprot

Q9UMR5

Synonyms:

Lysosomal thioesterase PPT2 (PPT-2) (EC 3.1.2) (S-thioesterase G14), PPT2

Immunogen:

Recombinant Human Lysosomal thioesterase PPT2 protein (28-302AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A