

TMEM37 Antibody, FITC conjugated



PACO64557

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Thought to stabilize the calcium channel in an inactivated (closed) state. Modulates calcium current when coexpressed with CACNA1G.

Gene ID:

TMEM37

Uniprot

Q8WXS4

Synonyms:

Voltage-dependent calcium channel gamma-like subunit (Neuronal voltage-gated calcium channel gamma-like subunit) (Transmembrane protein 37), TMEM37, PR

Immunogen:

Peptide sequence from Human Voltage-dependent calcium channel gamma-like subunit protein (4-20AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A