

# VPS26A Antibody, Biotin conjugated



PACO64562

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:****Gene ID:**

VPS26A

**Uniprot**

O75436

**Synonyms:**

Vacuolar protein sorting-associated protein 26A (Vesicle protein sorting 26A) (hVPS26), VPS26A, VPS26

**Immunogen:**

Recombinant Human Vacuolar protein sorting-associated protein 26A protein (21-180AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A