

VPS26A Antibody, Biotin conjugated



PACO64566

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

VPS26A

Uniprot

O75436

Synonyms:

Vacuolar protein sorting-associated protein 26A (Vesicle protein sorting 26A) (hVPS26), VPS26A, VPS26

Immunogen:

Recombinant Human Vacuolar protein sorting-associated protein 26A protein (21-180AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A