

CYR61 Antibody



PACO64575

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000, ELISA:1:10000

Protein Background:**Gene ID:**

CCN1

Uniprot

O00622

Synonyms:

3CH61 antibody; CCN 1 antibody; CCN family member 1 antibody; CCN1 antibody; CYR 61 antibody; Cyr61 antibody; CYR61 protein antibody; CYR61_HUMAN antibody;

Immunogen:

CYR61 Polyclonal Antibody detects endogenous levels of protein.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A