

CD20 Recombinant Antibody



RACO0031

Product Information

Size:

50ul

Reactivity:

Human

Source:

Human

Isotype:

Rabbit IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

This protein may be involved in the regulation of B-cell activation and proliferation.

Gene ID:

MS4A1

Uniprot

P11836

Synonyms:

B-lymphocyte antigen CD20, B-lymphocyte surface antigen B1, Bp35, Leukocyte surface antigen Leu-16, Membrane-spanning 4-domains subfamily A member 1, CD20, MS4A1, CD20

Immunogen:

A synthesized peptide.

Storage:

Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images

N/A

N/A