

# Histone H1.0 Recombinant Antibody



RACO0040

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Human

**Isotype:**

Rabbit IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Histones H1 are necessary for the condensation of nucleosome chains into higher-order structures. The H1F0 histones are found in cells that are in terminal stages of differentiation or that have low rates of cell division.

**Gene ID:**

H1F0

**Uniprot**

P07305

**Synonyms:**

Histone H1.0, Histone H1', Histone H1(0), Histone H1.0, N-terminally processed, H1F0, H1FV

**Immunogen:**

A synthesized peptide.

**Storage:**

Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---

N/A

N/A