

Phospho-PXN (Y118) Recombinant Antibody



RACO0108

Product Information

Size:

50ul

Reactivity:

Human

Source:

Human

Isotype:

Rabbit IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:5000

Protein Background:

Cytoskeletal protein involved in actin-membrane attachment at sites of cell adhesion to the extracellular matrix (focal adhesion).

Gene ID:

PXN

Uniprot

P49023

Synonyms:

Paxillin, PXN


Immunogen:

A synthesized peptide derived from human Phospho-PXN (Y118).

Storage:

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western Blot

Positive WB detected in (Hela whole cell lysate) A549 whole cell lysate)

HepG2 whole cell lysate) (treated with Calyculin A or EGF)

All lanes: Phospho-PXN antibody at 1: (3µg/ml)

Secondary

Goat polyclonal to rabbit IgG at 1:50000 dilution

Predicted band size: 68 KDa

Observed band size: 68 KDa