

# AXIN2 Recombinant Antibody



RACO0134

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Human

**Isotype:**

Rabbit IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Inhibitor of the Wnt signaling pathway. Down-regulates beta-catenin. Probably facilitate the phosphorylation of beta-catenin and APC by GSK3B (By similarity).

**Gene ID:**

AXIN2

**Uniprot**

Q9Y2T1

**Synonyms:**

Axin-2, Axin-like protein, Axil, Axis inhibition protein 2, Conductin, AXIN2

**Immunogen:**

A synthesized peptide derived from human AXIN2.

**Storage:**

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---

N/A

N/A