

# TBR1 Recombinant Antibody



RACO0150

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Human

**Isotype:**

Rabbit IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Probable transcriptional regulator involved in developmental processes. Required for normal brain development.

**Gene ID:**

TBR1

**Uniprot**

Q16650

**Synonyms:**

T-box brain protein 1, T-brain-1, TBR-1, TES-56, TBR1

**Immunogen:**

A synthesized peptide derived from human TBR1.

**Storage:**

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---

N/A

N/A