

PXN Recombinant Antibody



RACO0162

Product Information

Size:

50ul

Reactivity:

Human

Source:

Human

Isotype:

Rabbit IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Cytoskeletal protein involved in actin-membrane attachment at sites of cell adhesion to the extracellular matrix (focal adhesion).

Gene ID:

PXN

Uniprot

P49023

Synonyms:

Paxillin, PXN

Immunogen:

A synthesized peptide derived from human PXN.

Storage:

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images

N/A

N/A