

# IL10 Recombinant Antibody



RACO0176

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Human

**Isotype:**

Rabbit IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Inhibits the synthesis of a number of cytokines, including IFN-gamma, IL-2, IL-3, TNF and GM-CSF produced by activated macrophages and by helper T-cells.

**Gene ID:**

IL10

**Uniprot**

P22301

**Synonyms:**

Interleukin-10, IL-10, Cytokine synthesis inhibitory factor, CSIF, IL10

**Immunogen:**

A synthesized peptide derived from human IL10.

**Storage:**

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---

N/A

N/A