

GPX1 Recombinant Antibody



RACO0196

Product Information

Size:

50ul

Reactivity:

Human

Source:

Human

Isotype:

Rabbit IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Protects the hemoglobin in erythrocytes from oxidative breakdown.

Gene ID:

GPX1

Uniprot

P07203

Synonyms:

Glutathione peroxidase 1, GPx-1, GSHPx-1, Cellular glutathione peroxidase, GPX1

Immunogen:

A synthesized peptide derived from human GPX1.

Storage:

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images

N/A

N/A