

# CDC45 Recombinant Antibody



RACO0208

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Human

**Isotype:**

Rabbit IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Required for initiation of chromosomal DNA replication.

**Gene ID:**

CDC45

**Uniprot**

O75419

**Synonyms:**

Cell division control protein 45 homolog, PORC-PI-1, CDC45, CDC45L, CDC45L2, UNQ374/PRO710

**Immunogen:**

A synthesized peptide derived from human CDC45.

**Storage:**

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---

N/A

N/A