

# TGM2 Recombinant Antibody



RACO0343

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Homo sapiens (Human)

**Isotype:**

Rabbit IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:5000

**Protein Background:**

Catalyzes the cross-linking of proteins and the conjugation of polyamines to proteins.

**Gene ID:**

TGM2

**Uniprot**

P21980

**Synonyms:**

Protein-glutamine gamma-glutamyltransferase 2 (EC 2.3.2.13) (Tissue transglutaminase) (Transglutaminase C) (TG(C)) (TGC) (TGase C) (Transglutaminase H) (TGase H) (Transglutaminase-2) (TGase-2), TGM2

**Immunogen:**


A synthesized peptide derived from human TGM2.

**Storage:**

Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---



### Western Blot

Positive WB detected in( U87 whole cell lysate) All lanes: TGM2 antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1:50000 dilution

Predicted band size: 78, 62, 39 kDa

Observed band size: 78 kDa