

Technical Manual

ChromaDazzle Aconitase Activity Assay Kit

Catalogue Code: BA0070

Pack Size: 100 assays

Research Use Only

DESCRIPTION

ACONITASE (ACONITATE HYDRATASE) is an enzyme in the citric acid (TCA) cycle that catalyzes the conversion of citrate to isocitrate. The activity of aconitase depends largely upon the iron-sulfur $[\text{Fe}_4\text{S}_4]^{2+}$ cluster. Related diseases include aconitase deficiency (e.g. myopathy and exercise intolerance), Friedreich's ataxia and diabetes.

The Assay Genie ChromaDazzle Aconitase Activity Assay Kit measures the isocitrate generated as a product of the aconitase reaction. The isocitrate is then oxidized producing NADPH and the oxidation product. The NADPH converts the dye to an intense violet color with an absorption maximum at 565 nm. The increase in absorbance at 565 nm is directly proportional to aconitase activity.

KEY FEATURES

Fast and sensitive. Linear detection range (20 μL sample): 0.5 to 100 U/L for 20 min reaction.

Convenient and high-throughput. Homogeneous "mix-incubate-measure" type assay. Can be readily automated on HTS liquid handling systems for processing thousands of samples per day.

APPLICATIONS

Aconitase activity determination in biological samples (e.g. cell lysate, tissue homogenate, serum, etc.)

KIT CONTENTS (100 TESTS IN 96-WELL PLATES)

Assay Buffer:	10 mL	Enzyme A:	120 μL
NADP/MTT:	1 mL	Enzyme B:	120 μL
Standard:	1 mL	Substrate:	1 mL

Storage conditions. The kit is shipped on ice. Store all components at -20°C upon receiving. Shelf life: 6 months after receipt.

Precautions: reagents are for research use only. Normal precautions for laboratory reagents should be exercised while using the reagents. Please refer to Material Safety Data Sheet for detailed information.

PROCEDURES

Sample Preparation

Tissue: Prior to dissection, rinse tissue in phosphate buffered saline (pH 7.4) to remove blood. Homogenize tissue (50 mg) in $\sim 200 \mu\text{L}$ cold PBS. Centrifuge at $800 \times g$ for 10 min at 4°C . Remove supernatant for mitochondrial preparation.

Cell Lysate: Collect cells by centrifugation at $2,000 \times g$ for 5 min at 4°C . For adherent cells, do not harvest cells using proteolytic enzymes; rather use a rubber policeman. Homogenize or sonicate cells in an appropriate volume of cold PBS. Centrifuge at $800 \times g$ for 10 min at 4°C . Remove supernatant for mitochondrial preparation.

Mitochondrial Preparation: Centrifuge the removed supernatant at $20,000 \times g$ for 10 min at 4°C . Remove the supernatant and resuspend the pellet in cold PBS and sonicate for 20 seconds. The sample can be stored at -80°C for at least one month.

Reagent Preparation

Keep thawed Enzyme A and B on ice and equilibrate all other reagents to 25°C . Briefly centrifuge tubes before use.

Procedure using 96-well plate

1. **Standards.** Prepare 200 μL 5000 μM Premix by mixing 10 μL of the Standard (100 mM) and 190 μL distilled water. Dilute standards in 1.5-mL centrifuge tubes as described in the Table. Transfer 20 μL Standards into separate wells of a clear flat bottom 96-well plate.

No	Premix + H ₂ O	Isocitrate (μM)
1	100 μL + 0 μL	5000
2	60 μL + 40 μL	3000
3	30 μL + 70 μL	1500
4	0 μL + 100 μL	0

2. Transfer 20 μL of each sample into separate wells.

If any samples have high dehydrogenase activity or isocitrate levels, a sample blank will be needed. Transfer 20 μL of the sample into another well.

3. Prepare enough Working Reagent (WR) for sample and standard wells by mixing, for each well: 8 μL NADP/MTT Solution, 1 μL Enzyme A, 1 μL Enzyme B, 5 μL Substrate and 70 μL Assay Buffer. Fresh reconstitution of the WR is recommended.

If Sample Blanks are needed, prepare enough Blank Working Reagent (BWR) for the sample blank wells by mixing, for each well: 8 μL NADP/MTT Solution, 1 μL Enzyme B, and 75 μL Assay Buffer.

4. Add 80 μL WR or BWR to the corresponding wells. Tap plate briefly to mix.

5. Incubate at room temperature. Use a plate reader to read OD_{565nm} at 10 minutes and 30 minutes.

CALCULATION

Subtract blank value (water, #4) from the standard values at 10 min and plot the ΔOD against standard concentrations. Determine the slope and calculate the aconitase activity of the sample as follows

$$\text{Aconitase Activity} = \frac{\text{OD}_{30} - \text{OD}_{10}}{t \text{ (min)} \cdot \text{Slope } (\mu\text{M}^{-1})} \times n \quad (\text{U/L})$$


where OD₃₀, OD₁₀ are optical density values of the Sample at 30 and 10 minutes, respectively. *t* is the reaction time (20 minutes), Slope is the slope of the isocitrate standard curve and *n* is the sample dilution factor.

Unit Definition: 1 Unit (IU) of aconitase will catalyze the conversion of 1 μmole of citrate to isocitrate per minute at pH 7.4.

Note: if the calculated activity is higher than 100 U/L, dilute sample in water and repeat assay. Multiply the result by the dilution factor.

MATERIALS REQUIRED, BUT NOT PROVIDED

Pipetting devices and accessories (e.g. multi-channel pipettor), clear flat-bottom 96-well plates (e.g. VWR cat# 82050-760), centrifuge tubes and plate reader.



LITERATURE

1. Kennedy, MC et al. (1963) The Role of Iron in the Activation-Inactivation of Aconitase. J Biol Chem. 258(18):11098-105
2. Kennedy, MC et al. (1992) Purification and characterization of cytosolic aconitase from beef liver and its relationship to the iron responsive element binding protein. Proc Natl Acad Sci USA. 89(24):11730-34
3. Villafranca, JJ et al. (1974) The Mechanism of Aconitase Action. J Biol Chem. 249(19):6149-55.

Contact Details

Dublin, Ireland

Email: info@assaygenie.com

Web: www.assaygenie.com

Technical Support: Techsupport@assaygenie.com