



Technical Manual

Carbonyl Assay Kit

- **Catalogue Code: BA0206**
- **Size: 100 Assays**
- **Research Use Only**

Introduction

CARBONYL groups, such as ketones and aldehydes, are a common indicator of protein oxidation. Protein oxidation is caused as a result of exposure to reactive oxygen species (ROS) and is a common marker in various diseases as well as aging.

Assay Genies' Carbonyl Assay Kit is based on an improved method, where 2,4-dinitrophenylhydrazine (DNPH) reacts with carbonyl groups to produce a colored compound at 375 nm. The intensity of this colored compound is directly proportional to the carbonyl groups in the sample.

Key Features:

- **Fast and sensitive.** Use of 100 μ L sample. Linear detection range from 12 to 250 μ M carbonyl in 96-well plate assay.
- **Convenient.** The procedure involves adding a single working reagent, and reading the absorbance immediately after 30 min.
- **High-throughput.** Homogeneous "mix-incubate-measure" type assay.
Can be readily automated on HTS liquid handling systems for processing thousands of samples per day.

Application:

Direct Assays: Quantification of carbonyl groups (e.g. ketones, aldehydes) or protein carbonyls in biological samples (e.g. oxidized BSA, etc.)

Kit Contents:

Reagent: 12 mL

Standard: 50 μ L ((50 mM Carbonyl)

Additional materials needed:

- Pipetting devices
- centrifuge tubes
- clear flat-bottom 96-well plates (e.g. VWR cat# 82050-760)
- plate reader

Shipping Conditions:

The kit is shipped at RT. Store all components at -20°C upon receiving. Shelf life: 12 months after receipt.

Precautions:

Reagents are for research use only. Briefly centrifuge Standard tube before opening. Equilibrate all components to room temperature prior assay. Normal precautions for laboratory reagents should be exercised while using the reagents. Please refer to Material Safety Data Sheet for detailed information.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | **Web:** www.assaygenie.com

Copyright © 2020 Reagent Genie, All Rights Reserved. All information / detail is correct at time of going to print.

Assay Procedure

1. Standards. Prepare 1 mL of 250 μM Standard by mixing 5 μL of the 50 mM Standard and 995 μL dH₂O. Dilute standards in 1.5-mL centrifuge tubes as described in the Table.

	250 μM Standard + dH ₂ O	Carbonyl (μM)
1	500 μL + 0 μL	250
2	300 μL + 200 μL	150
3	150 μL + 350 μL	75
4	0 μL + 500 μL	0

2. Transfer 100 μL standards into separate wells of a clear, flat-bottom 96-well plate. Transfer 100 μL of each sample into separate wells.
3. Add 100 μL Reagent to the four Standards and the Sample Wells. Tap plate to mix briefly and thoroughly.
4. Incubate 30 min and read optical density at 375 nm (365-385 nm).

Calculations

Subtract the blank value (#4) from the standard values and plot the ΔOD against standard concentrations. Determine the slope and calculate the carbonyl concentration of Sample as follows:

$$[\text{Carbonyl}] = \frac{OD_{\text{Samples}} - OD_{\text{blank}}}{\text{Slope}(\mu\text{M}-1)} \times n (\mu\text{M})$$

Where

OD_{SAMPLE} and OD_{BLANK} are optical density readings of the Sample and Blank (Standard #4).
 n is the sample dilution factor.


Conversions: 10 μM Carbonyl equals 10 nmol/mL.

Note: it is recommended to normalize the carbonyl content to protein if measuring carbonyl protein.


Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | **Web:** www.assaygenie.com

Copyright © 2020 Reagent Genie, All Rights Reserved. All information / detail is correct at time of going to print.



Carbonyl Standard Curve



BSA Oxidation Time Course

BSA (1 mg/mL) was incubated with 12 mM acetaldehyde at 37°C for 1 week.

Carbonyl and protein content were measured. Carbonyl normalized to protein content.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | **Web:** www.assaygenie.com

Copyright © 2020 Reagent Genie, All Rights Reserved. All information / detail is correct at time of going to print.