

ARFGAP1 Antibody



PACO04717

Product Information

Size:

50ug

Gene ID:

ARFGAP1

Reactivity:

Human, Mouse, Rat, Monkey

Uniprot:

Q8N6T3

Source:

Rabbit

Synonyms:

ARFGAP1; ARF1GAP; ADP-ribosylation factor GTPase-activating protein 1; ARF GAP 1; ADP-ribosylation factor 1 GTPase-activating protein; ARF1 GAP; ARF1-directed GTPase-activating protein

Isotype:

IgG

Immunogen:

Synthesized peptide derived from the Internal region of Human ARF GAP1.

Applications:

WB, IHC, IF, ELISA;

Immunogen species:

Homo sapiens (Human)

Recommended dilutions:

WB: 1:500-1:2000

Storage Buffer:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

IHC: 1:100-1:300

IF: 1:200-1:1000

Storage:

Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.

ELISA: 1:20000

Form:

Liquid

Purification Method:

The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.


Conjugate:

Non-conjugated

Product Type:

Polyclonal Antibody

Product Images



Immunohistochemical analysis of paraffin-embedded human-liver, antibody was diluted at 1:200