



Recombinant Protein Technical Manual

Recombinant Human SPESP1 Protein (His Tag)

RPES0141

Product Data:

Product SKU: RPES0141

Size: 10µg

Species: Human

Expression host: Human Cells

Uniprot: Q6UW49

Protein Information:

Molecular Mass: 37.9 kDa

AP Molecular Mass: 58 kDa

Tag: C-6His

Bio-activity:

Purity: > 90 % as determined by reducing SDS-PAGE.

Endotoxin: < 1.0 EU per µg as determined by the LAL method.

Storage: Lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Shipping: This product is provided as lyophilized powder which is shipped with ice packs.

Formulation: Lyophilized from a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.

Reconstitution: Please refer to the printed manual for detailed information.

Application:

Synonyms: Sperm Equatorial Segment Protein 1; ESP; Equatorial Segment Protein; SP-ESP; Glycosylated 38 kDa Sperm Protein C-7/8; SPESP1

Immunogen Information:

Sequence: Tyr20-Tyr350

Background:

Sperm Equatorial Segment Protein 1 (SPESP1) is a member of the SPESP1 family. SPESP1 is highly expressed in the testis, where it is localized to the acrosome of postmeiotic stages of spermiogenesis; it is expressed at lower levels in the placenta and fetal lung. SPESP1 is involved in the multicellular organismal development. Disruption of SPESP1 leads to abnormal distribution of sperm proteins resulting in a detached membrane from the equatorial segment and less fertile sperm. SPESP1 may interact with IZUMO1 and MN9 antigen and it contains an N-glycosylation site as well as several cAMP-dependent kinase, protein kinase C, and casein kinase II consensus phosphorylation sites.