



Recombinant Protein Technical Manual

Recombinant Human SCF/c-Kit Ligand Protein (His Tag)(Active)

RPES0241

Product Data:

Product SKU: RPES0241

Size: 10µg

Species: Human

Expression host: Human Cells

Uniprot: P21583

Protein Information:

Molecular Mass: 22.0 kDa

AP Molecular Mass: 35 kDa

Tag: C-6His

Bio-activity: Measured in a cell proliferation assay using TF-1 human erythroleukemic cells. The ED50 for this effect is 20 ng/ml.

Purity: > 95 % as determined by reducing SDS-PAGE.

Endotoxin: < 1.0 EU per µg as determined by the LAL method.

Storage: Lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Shipping: This product is provided as lyophilized powder which is shipped with ice packs.

Formulation: Lyophilized from a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.

Reconstitution: Please refer to the printed manual for detailed information.

Application: Cell Culture

Synonyms: Kit Ligand; Mast Cell Growth Factor; MGF; Stem Cell Factor; SCF; c-Kit ligand; KITLG; MGF; SCF

Immunogen Information:

Sequence: Glu26-His214

Background:

Stem Cell Factor (SCF) is a hematopoietic growth factor that exerts its activity at the early stages of hematopoiesis. SCF stimulates the proliferation of myeloid, erythroid, and lymphoid progenitors in bone marrow cultures and has been shown to act synergistically with colony stimulating factors.