



Recombinant Protein Technical Manual

Recombinant Mouse DR6/TNFRSF21 Protein (His Tag)

RPES0323

Product Data:

Product SKU: RPES0323

Size: 10µg

Species: Mouse

Expression host: Human Cells

Uniprot: NP_055267.1

Protein Information:

Molecular Mass: 64.7 kDa

AP Molecular Mass:

Tag: C-6His

Bio-activity:

Purity: > 95 % as determined by SDS-PAGE

Endotoxin: < 1.0 EU per µg as determined by the LAL method.

Storage: Lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Shipping: This product is provided as lyophilized powder which is shipped with ice packs.

Formulation: Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.

Reconstitution: Please refer to it for detailed information.

Application:

Synonyms: Tumor necrosis factor receptor superfamily member 21;Death receptor 6;Tnfrsf21;CD358;BM-018;DR6

Immunogen Information:

Sequence: Gln42-His349

Background:

Tumor necrosis factor receptor superfamily member 21(DR6) is a single-pass type I membrane protein and contains 1 death domain and 4 TNFR-Cys repeats. The protein may activate NF-kappa-B and promote apoptosis and it may activate JNK and be involved in T-cell differentiation. It is required for both normal cell body death and axonal pruning. Trophic-factor deprivation triggers the cleavage of surface APP by beta-secretase to release sAPP-beta which is further cleaved to release an N-terminal fragment of APP (N-APP). N-APP binds TNFRSF21 triggering caspase activation and degeneration of both neuronal cell bodies (via caspase-3) and axons (via caspase-6).