

**Product Data:****Product SKU:** RPES0589**Size:** 50µg**Species:** Mouse**Expression host:** Baculovirus-Insect Cells**Uniprot:** P09411**Protein Information:****Molecular Mass:** 46 kDa**AP Molecular Mass:** 43 kDa**Tag:** C-His**Bio-activity:****Purity:** > 95 % as determined by SDS-PAGE**Endotoxin:** < 1.0 EU per µg of the protein as determined by the LAL method.**Storage:** Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles.**Shipping:** This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < -20°C.**Formulation:** Supplied as sterile 20mM Tris, 500mM NaCl, pH 7.4**Reconstitution:** Please refer to the printed manual for detailed information.**Application:****Synonyms:** Pgk

## Immunogen Information:

**Sequence:** Met 1-Val 417

## Background: