



Recombinant Protein Technical Manual

Recombinant Mouse Host Cell Factor 2/HCFC2 Protein (His Tag)

RPES0634

Product Data:

Product SKU: RPES0634

Size: 10µg

Species:

Expression host: E. coli

Uniprot: Q9D968

Protein Information:

Molecular Mass: 55.2 kDa

AP Molecular Mass: 55 kDa

Tag: N-6His

Bio-activity:

Purity: > 95 % as determined by reducing SDS-PAGE.

Endotoxin: < 1.0 EU per µg as determined by the LAL method.

Storage: Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles.

Shipping: This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < -20°C.

Formulation: Supplied as a 0.2 µm filtered solution of PBS, 2M Urea, pH 7.4.

Reconstitution: Please refer to the printed manual for detailed information.

Application:

Synonyms: Host cell factor 2 ;HCF-2;C2 factor;Hcfc2

Immunogen Information:

Sequence: Met 1-Glu497

Background: