



Recombinant Protein Technical Manual

Recombinant Mouse S100A4 Protein (His Tag)

RPES0759

Product Data:

Product SKU: RPES0759

Size: 10µg

Species: Mouse

Expression host: E. coli

Uniprot: P07091

Protein Information:

Molecular Mass: 12.5 kDa

AP Molecular Mass: 13 kDa

Tag: C-His

Bio-activity:

Purity: > 95% as determined by reducing SDS-PAGE.

Endotoxin: < 1.0 EU per µg as determined by the LAL method.

Storage: Lyophilized protein should be stored at < -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Shipping: This product is provided as lyophilized powder which is shipped with ice packs.

Formulation: Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.

Reconstitution: Please refer to the printed manual for detailed information.

Application:

Synonyms: Protein S100-A4; Metastasin; Metastatic cell protein; PEL98; Placental calcium-binding protein; Protein 18A2; Protein Mts1; S100 calcium-binding protein A4; S100a4; Capl; Mts1

Immunogen Information:

Sequence: Met1-Lys101

Background:

S100A4 is a member of the S100 family of proteins. The S100 family is further classified as a member of the EF-hand superfamily of Ca⁺⁺-binding proteins. These participate in both calcium-dependent and calcium-independent protein-protein interactions. The hallmark of this superfamily is the EF-hand motif that consists of a Ca⁺⁺-binding site flanked by two α -helices (helix E and helix F) that were originally identified in a right-handed model of carp muscle calcium-binding protein. Mouse S100A4 is 101 amino acids (aa) in length. It contains two EF hand domains, one between aa 12-47, and a second between aa 50-85. S100A4 activity has been associated with cell transformation. It seems likely this is either coincidental, or a consequence, rather than a cause of transformation.