



Recombinant Protein Technical Manual
Recombinant Mouse G-CSF/CSF3 Protein (His Tag)
RPES0938

Product Data:

Product SKU: RPES0938

Size: 10µg

Species: Mouse

Expression host: Human Cells

Uniprot: P09920

Protein Information:

Molecular Mass: 19.8 kDa

AP Molecular Mass: 25 kDa

Tag: C-6His

Bio-activity:

Purity: > 95 % as determined by SDS-PAGE

Endotoxin: < 1.0 EU per µg as determined by the LAL method.

Storage: Lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Shipping: This product is provided as lyophilized powder which is shipped with ice packs.

Formulation: Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4. .

Reconstitution: Please refer to the printed manual for detailed information.

Application:

Synonyms: Granulocyte colony-stimulating factor;Csf3;G-CSF

Immunogen Information:

Sequence: Val31-Ala208

Background:

Granulocyte colony-stimulating factor (G-CSF) is a growth factor and an essential cytokine which belongs to the IL-6 superfamily. Granulocyte/macrophage colony-stimulating factors are cytokines that act in hematopoiesis by controlling the production, differentiation, and function of 2 related white cell populations of the blood, the granulocytes and the monocytes-macrophages. G-CSF binding to its receptor G-CSF-R which belongs to the cytokine receptor type I family depends on the interaction of alpha-helical motifs of the former and two fibronectin type III as well as an immunoglobulin-like domain of the latter. G-CSF is a cytokine that have been demonstrated to improve cardiac function and perfusion in myocardial infarction. And it was initially evaluated as a stem cell mobilizer and erythropoietin as a cytoprotective agent.