



Recombinant Protein Technical Manual
Recombinant Human WWP2 Protein (His & GST Tag)
RPES1124

Product Data:

Product SKU: RPES1124

Size: 20µg

Species: Human

Expression host: Baculovirus-Insect Cells

Uniprot: O00308

Protein Information:

Molecular Mass: 126.7 kDa

AP Molecular Mass: 120 kDa

Tag: N-His & GST

Bio-activity:

Purity: > 90 % as determined by reducing SDS-PAGE.

Endotoxin: < 1.0 EU per µg as determined by the LAL method.

Storage: Lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Shipping: This product is provided as lyophilized powder which is shipped with ice packs.

Formulation: Lyophilized from sterile 20mM Tris, 500mM NaCl, pH 8.0, 10% glycerol

Reconstitution: Please refer to the printed manual for detailed information.

Application:

Synonyms: AIP2;WWp2-like

Immunogen Information:

Sequence: Met 1-Glu870

Background:

WWP2 contains 1 C2 domain, 1 HECT (E6AP-type E3 ubiquitin-protein ligase) domain and 4 WW domains. It is an E3 ubiquitin-protein ligase which accepts ubiquitin from an E2 ubiquitin-conjugating enzyme in the form of a thioester and then directly transfers the ubiquitin to targeted substrates. WWP2 can be detected in heart, throughout the brain, placenta, lung, liver, muscle, kidney and pancreas. It is also expressed in spleen and peripheral blood leukocytes. WWP2 polyubiquitinates POU5F1 by 'Lys-63'-linked conjugation and promotes it to proteasomal degradation; in embryonic stem cells (ESCs) the ubiquitination is proposed to regulate POU5F1 protein level. WWP2 ubiquitinates EGR2 and promotes it to proteasomal degradation; in T-cells the ubiquitination inhibits activation-induced cell death. It also ubiquitinates SLC11A2; the ubiquitination is enhanced by presence of NDFIP1 and NDFIP2. WWP2 ubiquitinates RPB1 and promotes it to proteasomal degradation.