



Recombinant Protein Technical Manual  
Recombinant Mouse CHODL/Chondrolectin Protein  
(Fc Tag)  
RPES1982

#### Product Data:

**Product SKU:** RPES1982

**Size:** 20µg

**Species:** Mouse

**Expression host:** HEK293 Cells

**Uniprot:** EDK98299.1

#### Protein Information:

**Molecular Mass:** 48.9 kDa

**AP Molecular Mass:** 60 kDa

**Tag:** C-Fc

**Bio-activity:**

**Purity:** > 95 % as determined by SDS-PAGE

**Endotoxin:** < 1.0 EU per µg of the protein as determined by the LAL method.

**Storage:** Lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

**Shipping:** This product is provided as lyophilized powder which is shipped with ice packs.

**Formulation:** Lyophilized from sterile PBS, pH 7.4

**Reconstitution:** Please refer to the printed manual for detailed information.

**Application:**

**Synonyms:** 3110074E07Rik;MT75;PRED12

## Immunogen Information:

**Sequence:** Met 1-Asn 216

## Background: