



Recombinant Protein Technical Manual
Recombinant Human DNMT2/TRDMT1 Protein (GST
Tag)
RPES2635

Product Data:

Product SKU: RPES2635

Size: 20µg

Species: Human

Expression host: Baculovirus-Insect Cells

Uniprot: NP_004403.1

Protein Information:

Molecular Mass: 71 kDa

AP Molecular Mass: 60 kDa

Tag: N-GST

Bio-activity:

Purity: > 94 % as determined by reducing SDS-PAGE.

Endotoxin: < 1.0 EU per µg as determined by the LAL method.

Storage: Lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Shipping: This product is provided as lyophilized powder which is shipped with ice packs.

Formulation: Lyophilized from sterile 50mM Tris, 100mM NaCl, 0.5mM GSH, 0.5mM PMSF, pH 8.0

Reconstitution: Please refer to the printed manual for detailed information.

Application:

Synonyms: DMNT2;DNMT2;MHSIIP;PUMET;RNMT1

Immunogen Information:

Sequence: Met 1-Glu 391

Background:

DNMT2, also known as tRNA (cytosine-5-)-methyltransferase, DNA methyltransferase homolog HsaIIP, and TRDMT1, is a member of the DNA methyltransferase family of enzymes. DNMT2 enzymes have been widely conserved during evolution and contain all of the signature motifs of DNA (cytosine-5)-methyltransferases. It contains all 10 sequence motifs that are conserved among m(5)C MTases, including the consensus S-adenosyl-L-methionine-binding motifs and the active site ProCys dipeptide, and its structure is very similar to prokaryotic DNA methyltransferases. DNMT2 has close homologs in plants, insects and *Schizosaccharomyces pombe*, but no related sequence can be found in the genomes of *Saccharomyces cerevisiae* or *Caenorhabditis elegans*. While the biological function of DNMT2 is not yet known, the strong binding to DNA suggests that DNMT2 may mark specific sequences in the genome by binding to DNA through the specific target-recognizing motif. However, the DNA methyltransferase activity of these proteins is comparatively weak and their biochemical and functional properties remain enigmatic. Recent evidence now shows that Dnmt2 has a novel tRNA methyltransferase activity, raising the possibility that the biological roles of these proteins might be broader than previously thought.