



Recombinant Protein Technical Manual
Recombinant Human Apolipoprotein E/ApoE Protein
(His Tag)
RPES2909

Product Data:

Product SKU: RPES2909

Size: 10µg

Species: Human

Expression host: Human Cells

Uniprot: P02649

Protein Information:

Molecular Mass: 35.3 kDa

AP Molecular Mass: 39 kDa

Tag: C-6His

Bio-activity:

Purity: > 95 % as determined by reducing SDS-PAGE.

Endotoxin: < 1.0 EU per µg as determined by the LAL method.

Storage: Lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Shipping: This product is provided as lyophilized powder which is shipped with ice packs.

Formulation: Lyophilized from a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.

Reconstitution: Please refer to the printed manual for detailed information.

Application:

Synonyms: Apolipoprotein E;APOE;Apo-E

Immunogen Information:

Sequence: Lys19-His317

Background:

ApoE, a glycoprotein, is a structural component of very low density lipoprotein (vLDL) synthesized by the liver and intestinally synthesized chylomicrons . ApoE is also a constituent of a subclass of high density of lipoproteins (HDL) involved in cholesterol transport . ApoE mediates high affinity binding of chylomicrons and vLDL particles to the LDL receptor, allowing for specific uptake of these particles by the liver, preventing the accumulation of cholesterol rich particles in the plasma . Apolipoprotein E combines with fats (lipids) in the body to form molecules called lipoproteins and Apolipoprotein E is a major component of a specific type of lipoprotein called very low-density lipoproteins (VLDLs).