



Recombinant Protein Technical Manual

**Recombinant Mouse BAFFR/TNFRSF13C/CD268
Protein (Fc Tag)(Active)**
RPES2951

Product Data:

Product SKU: RPES2951

Size: 50µg

Species: Mouse

Expression host: HEK293 Cells

Uniprot: Q9D8D0

Protein Information:

Molecular Mass: 33.7 kDa

AP Molecular Mass: 40-45 kDa

Tag: C-Fc

Bio-activity: Immobilized human BAFF at 10 µg/ml (100 µl/well) can bind mouse BAFFR-Fc, The EC50 of mouse BAFFR-Fc is 0.14-0.32 µg/ml.

Purity: > 90 % as determined by SDS-PAGE

Endotoxin: < 1.0 EU per µg of the protein as determined by the LAL method.

Storage: Lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Shipping: This product is provided as lyophilized powder which is shipped with ice packs.

Formulation: Lyophilized from sterile PBS, pH 7.4

Reconstitution: Please refer to the printed manual for detailed information.

Application: Functional ELISA

Synonyms: 2010006P15Rik;BAFF-R;Baffr;Bcmd;Bcmd;Bcmd1;Lvis22

Immunogen Information:

Sequence: Met 1-Ala 71

Background:

Tumor necrosis factor receptor superfamily, member 13C (TNFRSF13C) also known as B-cell-activating factor receptor (BAFFR) and CD268 antigen, is a member of the tumor necrosis factor receptor superfamily. A tumor necrosis factor receptor (TNFR), or death receptor, is a trimeric cytokine receptor that binds tumor necrosis factors (TNF). The receptor cooperates with an adaptor protein which is important in determining the outcome of the response. Members of the TNF receptor superfamily (TNFRSF) have crucial roles in both innate and adaptive immunity and in cellular apoptosis process. Apoptosis is a cell suicide mechanism that enables metazoans to control cell number in tissues and to eliminate individual cells that threaten the animal's survival. Certain cells have unique sensors, termed death receptors or tumour necrosis factor (TNFR), on their surface. Tumour necrosis factors (TNFR) detect the presence of extracellular death signals and, in response, they rapidly ignite the cell's intrinsic apoptosis machinery. It has been proposed that abnormally high levels of BAFFR/TNFRSF13C (CD268) may contribute to the pathogenesis of autoimmune diseases by enhancing the survival of autoreactive B cells.