



Recombinant Protein Technical Manual

Recombinant Mouse BTLA/CD272 Protein (Fc Tag)(Active)
RPES2999

Product Data:

Product SKU: RPES2999

Size: 50µg

Species: Mouse

Expression host: HEK293 Cells

Uniprot: Q7TSA3

Protein Information:

Molecular Mass: 437 kDa

AP Molecular Mass: 57-62 kDa

Tag: C-Fc

Bio-activity: Immobilized mouse BTLA-Fc at 10 µg/ml (100 µl/well) can bind biotinylated mouse HVEM-Fch, The EC50 of biotinylated mouse HVEM-Fch is 15.7-36.7 ng/ml.

Purity: > 90 % as determined by SDS-PAGE

Endotoxin: < 1.0 EU per µg of the protein as determined by the LAL method.

Storage: Lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Shipping: This product is provided as lyophilized powder which is shipped with ice packs.

Formulation: Lyophilized from sterile PBS, pH 7.4

Reconstitution: Please refer to the printed manual for detailed information.

Application: Functional ELISA

Synonyms: A630002H24

Immunogen Information:

Sequence: Met1-Gly176

Background:

BTLA is an inhibitory molecule which belongs to the Ig superfamily. It down-modulates immune responses. As such, reagents that regulate the binding of BTLA to its ligand or alter BTLA signaling have significant therapeutic promise. BTLA is crucial to understand the mechanism(s) of action of these antibodies before attempting clinical applications. BTLA is not expressed by naive T cells, but it is induced during activation and remains expressed on T helper type 1 (T(H)1) but not T(H)2 cells. BTLA is a third inhibitory receptor on T lymphocytes with similarities to cytotoxic T lymphocyte-associated antigen 4 (CTLA-4) and programmed death 1 (PD).