



# Recombinant Protein Technical Manual

## Recombinant Human FAM3D Protein (His Tag)

RPES3203

### Product Data:

**Product SKU:** RPES3203

**Size:** 10µg

**Species:** Human

**Expression host:** Human Cells

**Uniprot:** Q96BQ1

### Protein Information:

**Molecular Mass:** 23.1 kDa

**AP Molecular Mass:** 26 kDa

**Tag:** C-6His

**Bio-activity:**

**Purity:** > 95 % as determined by reducing SDS-PAGE.

**Endotoxin:** < 1.0 EU per µg as determined by the LAL method.

**Storage:** Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles.

**Shipping:** This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < -20°C.

**Formulation:** Supplied as a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.2.

**Reconstitution:** Please refer to the printed manual for detailed information.

**Application:**

**Synonyms:** Protein FAM3D; FAM3D

## Immunogen Information:

**Sequence:** Tyr26-Phe224

## Background:

Protein FAM3D is a novel cytokine-like protein that belongs to the FAM3 family. Human FAM3D is synthesized as a 224 amino acid precursor that contains a 25 amino acid signal sequence and a 199 amino acid mature chain. FAM3D is identified based on structural, but not sequence, homology to short chain cytokines including IL-2, IL-4 and GM-CSF. FAM3 proteins are four helix bundle cytokines with four conserved cysteines in all members (FAM3A-D). FAM3B is highly expressed in alpha and beta cells of the pancreas and is being investigated as a potential contributor to beta cell death and development of Type I Diabetes.