



Recombinant Protein Technical Manual
Recombinant Human Alpha-Galactosidase A/GLA
Protein (His Tag)
RPES3204

Product Data:

Product SKU: RPES3204

Size: 10µg

Species: Human

Expression host: Human Cells

Uniprot: P06280

Protein Information:

Molecular Mass: 46.4 kDa

AP Molecular Mass: 51 kDa

Tag: C-6His

Bio-activity:

Purity: > 95 % as determined by reducing SDS-PAGE.

Endotoxin: < 1.0 EU per µg as determined by the LAL method.

Storage: Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles.

Shipping: This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < -20°C.

Formulation: Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, pH 8.0.

Reconstitution: Please refer to the printed manual for detailed information.

Application:

Synonyms: Alpha-Galactosidase A; Alpha-D-Galactosidase A; Alpha-D-Galactoside Galactohydrolase; Melibiase; Agalsidase; GLA;GLAL

Immunogen Information:

Sequence: Leu32-Leu429

Background:

α -Galactosidase A is a homodimeric glycoprotein that belongs to the glycosyl hydrolase 27 family. It is a lysosomal enzyme and used as a long-term enzyme replacement therapy in patients with a confirmed diagnosis of Fabry disease. α -Galactosidase A can hydrolyze terminal α -galactosyl moieties from glycolipids and glycoproteins and catalyze the hydrolysis of melibiose into galactose and glucose. Defects α -Galactosidase A are the cause of Fabry disease (FD) which is a rare X-linked sphingolipidosis disease with glycolipid accumulates in many tissues. The disease consists of an inborn error of glycosphingolipid catabolism. FD patients show systemic accumulation of globotriaoslyceramide (Gb3) and related glycosphingolipids in the plasma and cellular lysosomes throughout the body. Patients may show ocular deposits, febrile episodes, and burning pain in the extremities. Death results from renal failure, cardiac or cerebral complications of hypertension or other vascular disease.