



Recombinant Protein Technical Manual

Recombinant Rat HAAO/3-HAO Protein (His Tag)

RPES3297

Product Data:

Product SKU: RPES3297

Size: 20µg

Species: Rat

Expression host: E. coli

Uniprot: NP_064461.1

Protein Information:

Molecular Mass: 34.8 kDa

AP Molecular Mass:

Tag: N-His

Bio-activity:

Purity: > 90 % as determined by SDS-PAGE

Endotoxin: Please contact us for more information.

Storage: Lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Shipping: This product is provided as lyophilized powder which is shipped with ice packs.

Formulation: Lyophilized from sterile PBS, pH 7.4

Reconstitution: Please refer to the printed manual for detailed information.

Application:

Synonyms: HAAO;3-HAO

Immunogen Information:

Sequence: Met1-Gly286

Background: