



Recombinant Protein Technical Manual

Recombinant Human LMIR2/CD300C Protein (Fc Tag)

RPES3526

Product Data:

Product SKU: RPES3526

Size: 10µg

Species: Human

Expression host: Human Cells

Uniprot: Q08708

Protein Information:

Molecular Mass: 44.9 kDa

AP Molecular Mass: 60-90 kDa

Tag: C-Fc

Bio-activity:

Purity: > 90 % as determined by reducing SDS-PAGE.

Endotoxin: < 1.0 EU per µg as determined by the LAL method.

Storage: Lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Shipping: This product is provided as lyophilized powder which is shipped with ice packs.

Formulation: Lyophilized from a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.

Reconstitution: Please refer to the printed manual for detailed information.

Application:

Synonyms: CMRF35-Like Molecule 6; CLM-6; CD300 Antigen-Like Family Member C; CMRF35-A1; CMRF-35; Immunoglobulin Superfamily Member 16; IgSF16; CD300c; CD300C; CMRF35; CMRF35A; CMRF35A1; IGSF16

Immunogen Information:

Sequence: Gly21-Arg183

Background:

CD300C is a single-pass type I membrane protein which belongs to the immunoregulatory signaling (IRS) family. CD300C contains one Ig-like V-type domain and is present on the surface of natural killer cells, granulocytes, most myeloid cells, dendritic cells, and a subpopulation of T and B lymphocytes. The CD300C (CMRF-35A) and CD300A (CMRF-35H) molecules are homologous leukocyte surface proteins. CD300a and CD300C play an important role in the cross-regulation of TNF-alpha and IFN-alpha secretion from pDCs. CD300A and CD300C are indistinguishable on the surface of NK cells. The ligand for CD300C is presently unknown.