



Recombinant Protein Technical Manual
Recombinant Human XRP2/RP2 Protein (GST Tag)
RPES3642

Product Data:

Product SKU: RPES3642

Size: 10µg

Species: Human

Expression host: Baculovirus-Insect Cells

Uniprot: NP_008846.2

Protein Information:

Molecular Mass: 66 kDa

AP Molecular Mass: 66 kDa

Tag: N-GST

Bio-activity:

Purity: > 95 % as determined by reducing SDS-PAGE.

Endotoxin: < 1.0 EU per µg as determined by the LAL method.

Storage: Lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Shipping: This product is provided as lyophilized powder which is shipped with ice packs.

Formulation: Lyophilized from sterile 50mM Tris, 100mM NaCl, 1mM GSH 0.5mM EDTA, 0.5mM PMSF pH 8.0

Reconstitution: Please refer to the printed manual for detailed information.

Application:

Synonyms: DELXp11.3;KIAA0215;NM23-H10;NME10;RP2;TBCCD2;XRP2

Immunogen Information:

Sequence: Met 1-Thr 350

Background:

XRP2, also known as Protein XRP2 and RP2, is a member of the TBCC (tubulin cofactor C) family and contains one C-CAP/cofactor C-like domain. This protein is encoded by the RP2 gene in humans. XRP2 stimulates the GTPase activity of tubulin, but does not enhance tubulin heterodimerization. XRP2 acts as guanine nucleotide dissociation inhibitor for ARL3. Defects in RP2 gene are the cause of retinitis pigmentosa type 2 (RP2), also known as X-linked retinitis pigmentosa 2 (XLRP-2). It leads to degeneration of retinal photoreceptor cells. Patients typically have night vision blindness and loss of midperipheral visual field. As their condition progresses, they lose their far peripheral visual field and eventually central vision as well.