



Recombinant Protein Technical Manual
Recombinant Human RAB31/Rab22B Protein (His
Tag)
RPES3727

Product Data:

Product SKU: RPES3727

Size: 20µg

Species: Human

Expression host: E. coli

Uniprot: NP_006859.2

Protein Information:

Molecular Mass: 23.9 kDa

AP Molecular Mass: 26 kDa

Tag: N-His

Bio-activity:

Purity: > 90 % as determined by reducing SDS-PAGE.

Endotoxin: Please contact us for more information.

Storage: Lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Shipping: This product is provided as lyophilized powder which is shipped with ice packs.

Formulation: Lyophilized from sterile 50 mM Tris, pH 8.0

Reconstitution: Please refer to the printed manual for detailed information.

Application:

Synonyms: Rab22B

Immunogen Information:

Sequence: Met 1-Cys195

Background:

Small GTP-binding proteins of the RAB family, such as RAB31, play essential roles in vesicle and granule targeting (Bao et al. , 2002 [PubMed 11784320])