



Recombinant Protein Technical Manual

Recombinant Human SMR3B Protein (His Tag)

RPES3951

Product Data:

Product SKU: RPES3951

Size: 20µg

Species: Human

Expression host: HEK293 Cells

Uniprot: NP_006676.1

Protein Information:

Molecular Mass: 7.3 kDa

AP Molecular Mass: 8 kDa

Tag: C-His

Bio-activity:

Purity: > 95 % as determined by reducing SDS-PAGE.

Endotoxin: < 1.0 EU per µg of the protein as determined by the LAL method.

Storage: Lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Shipping: This product is provided as lyophilized powder which is shipped with ice packs.

Formulation: Lyophilized from sterile PBS, pH 7.4

Reconstitution: Please refer to the printed manual for detailed information.

Application:

Synonyms: P-B;PBII;PRL3;PROL3;SMR1B;SMR3B

Immunogen Information:

Sequence: Met 1-Pro79

Background: