



Recombinant Protein Technical Manual

Recombinant Human PNLIP Protein (His Tag)

RPES4625

Product Data:

Product SKU: RPES4625

Size: 20µg

Species: Human

Expression host: HEK293 Cells

Uniprot: P16233

Protein Information:

Molecular Mass: 51 kDa

AP Molecular Mass: 51 kDa

Tag: C-His

Bio-activity:

Purity: > 95 % as determined by reducing SDS-PAGE.

Endotoxin: < 1.0 EU per µg as determined by the LAL method.

Storage: Lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Shipping: This product is provided as lyophilized powder which is shipped with ice packs.

Formulation: Lyophilized from sterile PBS, pH 7.4

Reconstitution: Please refer to the printed manual for detailed information.

Application:

Synonyms: PL;PNLIPD;PTL

Immunogen Information:

Sequence: Met 1-Cys465

Background:

PNLIP is an enzyme which belongs to the lipase family. Secreted from the pancreas, PNLIP is the primary lipase that hydrolyzes dietary fat molecules in the human digestive system, converting triglyceride substrates found in ingested oils to monoglycerides and free fatty acids. Bile salts secreted from the liver and stored in gallbladder are released into the duodenum where they coat and emulsify large fat droplets into smaller droplets, thus increasing the overall surface area of the fat, which allows the lipase to break apart the fat more effectively. The resulting monomers (2 free fatty acids and one 2-monoacylglycerol) are then moved by way of peristalsis along the small intestine to be absorbed into the lymphatic system by a specialized vessel called a lacteal.