



Recombinant Protein Technical Manual
Recombinant Human CD157/BST1 Protein (Fc Tag)
RPES4750

Product Data:

Product SKU: RPES4750

Size: 10µg

Species: Human

Expression host: Human Cells

Uniprot: Q10588

Protein Information:

Molecular Mass: 57.1 kDa

AP Molecular Mass: 75 kDa

Tag: C-Fc

Bio-activity:

Purity: > 95% as determined by reducing SDS-PAGE.

Endotoxin: < 1.0 EU per µg as determined by the LAL method.

Storage: Lyophilized protein should be stored at < -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Shipping: This product is provided as lyophilized powder which is shipped with ice packs.

Formulation: Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.

Reconstitution: Please refer to the printed manual for detailed information.

Application:

Synonyms: ADP-ribosyl cyclase 2; Cyclic ADP-ribose hydrolase 2; BST; Bone marrow stromal antigen 1; Bone marrow stromal cell antigen 1

Immunogen Information:

Sequence: Gly29-Lys292

Background:

CD38, also called ADP-ribosyl cyclase, is a Type II integral membrane protein with 301 amino acids in length that belongs to the ADP-ribosyl cyclase family. It synthesizes the second messengers cyclic ADP-ribose and nicotinate-adenine dinucleotide phosphate, the former a second messenger for glucose-induced insulin secretion. And also moonlights as a receptor in cells of the immune system. CD38 is expressed in B and T lymphocytes, osteoclasts, and in cardiac, pancreatic, liver and kidney cells. Through its production of cyclic ADP-ribose, CD38 modulates calcium-mediated signal transduction in many types of cells, including neutrophils and pancreatic beta cells.