



Recombinant Protein Technical Manual
Recombinant Mouse CD79A/Ig-alpha Protein (ECD,
Fc Tag)
RPES4940

Product Data:

Product SKU: RPES4940

Size: 20µg

Species: Mouse

Expression host: HEK293 Cells

Uniprot: NP_031681.2

Protein Information:

Molecular Mass: 39 kDa

AP Molecular Mass:

Tag: C-Fc

Bio-activity:

Purity: > 90 % as determined by SDS-PAGE

Endotoxin: < 1.0 EU per µg of the protein as determined by the LAL method.

Storage: Lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Shipping: This product is provided as lyophilized powder which is shipped with ice packs.

Formulation: Lyophilized from sterile PBS, pH 7.4

Reconstitution: Please refer to the printed manual for detailed information.

Application:

Synonyms: Ig-alpha;lga;lgalpha;Ly-54;Ly54;mb

Immunogen Information:

Sequence: Met1-Arg137

Background: