

Recombinant Protein Technical Manual

Recombinant Mouse Syndecan-4/SDC4 Protein (aa 146, His Tag) RPES5022

Product Data:

Product SKU: RPES5022

Size: 50µg

Species: Mouse

Expression host: HEK293 Cells

Uniprot: 035988

Dura ta ta	
Protein	KIT IONY

Molecular Mass:	14.9 kDa
AP Molecular Mass:	28 kDa
Tag:	C-His
Bio-activity:	
Purity:	> 95 % as determined by SDS-PAGE
Endotoxin:	< 1.0 EU per μg of the protein as determined by the LAL method.
Storage:	Lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping:	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation:	Lyophilized from sterile PBS, pH 7.4
Reconstitution:	Please refer to the printed manual for detailed information.
Application:	
Synonyms:	AA959608;AW108331;ryudocan;Synd4;syndecan-4;SDC4;SYND4;Ryudocan core protein

Sequence: Met 1-Val 146

Background:

SDC4 (Syndecan-4), also known as Syn4, is a transmembrane heparan sulfate proteoglycan that co-operates with integrins during cell-matrix interactions for the assembly of focal adhesions and actin stress fibers and in the phosphorylation of focal adhesion kinase (FAK) on Tyr397. Syndecan-4 plays roles in the formation of focal adhesions and stress fibers. The cytoplasmic domain of syndecan-4 interacts with a number of signalling and structural proteins, and both extracellular and cytoplasmic domains are necessary for regulated activation of associated transmembrane receptors. Syndecan-4/SDC4 is a heparan sulfate proteoglycan and works as a coreceptor for various growth factors. SDC4 deficiency limits neointimal formation after vascular injury by regulating vascular smooth muscle cells (VSMCs) proliferation and vascular progenitor cells (VPCs) mobilization. Therefore, SDC4 may be a novel therapeutic target for preventing arterial restenosis after angioplasty.