

# Human APOA1 Recombinant Protein



RPPB0047

---

## Product Information    Protein Information

**Product SKU:**

RPPB0047

**Accession:**

P02647

**Host:**

Human HDL.

**Protein description:**

APOA1 Human isolated from Human HDL is a single, glycosylated, polypeptide chain having a molecular mass of 28.3kDa. APOA1 is purified using delipidation and gel permeation chromatographic technique.

**Appearance:**

Sterile filtered colorless solution.

**Synonyms:**

Apolipoprotein A-I, Apo-AI, ApoA-I, APOA1, MGC117399.

**Formulation:**

The APOA1 1mg/ml solution contains 10mM Ammonium Bicarbonate at pH7.4.

**Purity:**

Greater than 95.0% as determined by SDS-PAGE.

**Stability:**

APOA1 although stable at 4°C for 1 week, should be stored below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent freeze thaw cycles.