

Human APO D Recombinant Protein



RPPB0052

Product Information Protein Information

Product SKU:

RPPB0052

Accession:

P05090

Host:

Escherichia Coli.

Protein description:

APO-D Human Recombinant full length protein expressed in E.coli, shows a 48 kDa band on SDS-PAGE. The APO-D is purified by proprietary chromatographic techniques.

Appearance:

Sterile Filtered clear solution.

Synonyms:

Apolipoprotein D, Apo-D, ApoD.

Formulation:

APO-D protein at 100µg/ml in 50mM Tris-Acetate, pH7.5, 1mM EDTA and 20% Glycerol.

Stability:

Store vial at -20°C to -80°C. When stored at the recommended temperature, this protein is stable for 12 months.