

# Human Clusterin Recombinant Protein



RPPB0060

---

## Product Information    Protein Information

**Product SKU:**

RPPB0060

**Appearance:**

Filtered White lyophilized (freeze-dried) powder.

**Accession:**

P10909

**Synonyms:**

CLI, AAG4, KUB1, SGP2, SGP-2, SP-40, TRPM2, MGC24903, Clusterin, Apolipoprotein J, Apo-J.

**Host:**

Plasma.

**Formulation:**

Human native Clusterin was filtered (0.4µm) and lyophilized from 0.5mg/ml in 0.1M phosphate buffer, 0.15M NaCl pH 7.5.

**Purity:**

Greater than 95% as determined by SDS PAGE.

**Solubility:**

It is recommended to add deionized H<sub>2</sub>O to prepare a working stock solution of 0.5mg/ml and let the lyophilized pellet dissolve completely. Product is not sterile! Please filter the product by an appropriate sterile filter before using it in the cell culture.

**Stability:**

Store lyophilized Clusterin at -20°C. Aliquot the product after reconstitution to avoid repeated freezing/thawing cycles. Reconstituted protein can be stored at 4°C for a limited period of time; it does not show any change after two weeks at 4°C.