

Human Clusterin Recombinant Protein



RPPB0061

Product Information Protein Information

Product SKU:

RPPB0061

Appearance:

Filtered White lyophilized (freeze-dried) powder.

Accession:

P10909

Synonyms:

CLI, AAG4, KUB1, SGP2, SGP-2, SP-40, TRPM2, MGC24903, Clusterin, Apolipoprotein J, Apo-J.

Host:

Human Serum.

Formulation:

Human native Clusterin was filtered (0.4 μ m) and lyophilized from 0.5mg/ml in phosphate buffer saline, pH 7.5 and 5% (w/v) trehalose.

Purity:

Greater than 90% as determined by SDS PAGE.

Solubility:

It is recommended to add deionized H₂O to prepare a working stock solution of 0.5mg/ml and let the lyophilized pellet dissolve completely. Product is not sterile! Please filter the product by an appropriate sterile filter before using it in the cell culture.

Stability:

Store lyophilized Clusterin at -20°C. Aliquot the product after reconstitution to avoid repeated freezing/thawing cycles. Reconstituted protein can be stored at 4°C for a limited period of time; it does not show any change after two weeks at 4°C.