

# Mouse FGF 15 Recombinant Protein



RPPB0228

---

## Product Information    Protein Information

**Product SKU:**

RPPB0228

**Accession:**

O35622

**Host:**

E. coli

**Protein description:**

FGF-15 Protein is 23.78 kDa protein containing 203 amino acid residues including a 10 a.a. N-Terminal His-tag.

**Appearance:**

Filtered colorless solution.

**Synonyms:**

Fibroblast Growth Factor 15, FGF19.

**Formulation:**

The FGF-15 filtered (0.4 $\mu$ m) solution in 20mM TRIS and 50mM NaCl, pH 8.5.

**Purity:**

Greater than 90.0% as determined by SDS-PAGE.

**Stability:**

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

**Amino Acid Sequence:**

MKHHHHHHAS RPLAQQSQSV SDEDPLFLYG WGKITRLQYL YSAGPYVSNC FLIRSDGSV DCEEDQNERN  
LLEFRAVALK TIAIKDVSSV RYLCMSADGK IYGLIRYSEE DCTFREEMDC LGYNQYRSMK HHLHIIIFIQA  
KPREQLQDQK PSNFIPVFHR SFFETGDQLR SKMFSLPLES DSMDFRMVE DVDHLVKSPS FQK.