

# Human IGF1 Recombinant Protein



RPPB0423

---

## Product Information    Protein Information

**Product SKU:**

RPPB0423

**Accession:**

P05019

**Host:**

Escherichia Coli.

**Protein description:**

IGF1 A70T Human Recombinant produced in E.Coli is a single, non-glycosylated polypeptide chain containing 70 amino acids and having a molecular mass of Approximately 7.7kDa. The IGF1 A70T is purified by proprietary chromatographic techniques.

**Appearance:**

Sterile Filtered White lyophilized (freeze-dried) powder.

**Synonyms:**

Somatomedin C, IGF-I, IGFI, IGF1, IGF-IA, Mechano growth factor, MGF.

**Formulation:**

IGF1 A70T Lyophilized from a 0.2 µm filtered concentrated solution in 20 mM PB and 150 mM NaCl, pH 6.0.

**Purity:**

Greater than 95.0% as determined by SDS-PAGE.

**Solubility:**

It is recommended to reconstitute the lyophilized IGF1 A70T in sterile 18MΩ-cm H<sub>2</sub>O not less than 100µg/ml, which can then be further diluted to other aqueous solutions.

**Stability:**

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.

**Amino Acid Sequence:**

GPETLCGAEL VDALQFVCGD RGFYFNKPTG YGSSRRAPQ TGIVDECCFR SCDLRRLEMY CAPLPAKST.