

Human IL 1RA His Recombinant Protein



RPPB0466

Product Information Protein Information

Product SKU:

RPPB0466

Accession:

P18510

Host:

Escherichia Coli.

Protein description:

IL1Ra Human Recombinant produced in E.Coli is single, a non-glycosylated, Polypeptide chain containing 152 amino acids fragment (26-177) and having a total molecular mass of 21.63kDa with an amino-terminal hexahistidine tag. The IL-1Ra His is purified by proprietary chromatographic techniques.

Appearance:

Sterile Filtered clear solution.

Synonyms:

IRAP, IL1F3, IL1RA; IL-1ra3, ICIL-1RA, IL1RN, IL1 inhibitor, IL-1ra, MGC10430.

Formulation:

Interleukin-1Ra His-tag protein is supplied in 20mM Tris HCl (pH 8) and 50% glycerol.

Purity:

Greater than 95.0% as determined by SDS-PAGE.

Stability:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Please avoid freeze thaw cycles.