

# Human IL 7 Recombinant Protein



RPPB0540

---

## Product Information    Protein Information

**Product SKU:**

RPPB0540

**Accession:**

P13232

**Host:**

Escherichia Coli.

**Protein description:**

IL-7 Human Recombinant produced in E.Coli is single, a non-glycosylated, Polypeptide chain containing 152 amino acids fragment (26-177) and having a total molecular mass of 21.97 kDa with an amino-terminal hexahistidine tag. The IL-7 His-Tag protein is purified by proprietary chromatographic techniques.

**Appearance:**

Sterile Filtered clear solution.

**Synonyms:**

Lymphopoietin 1 (LP-1), pre-B cell factor, IL-7.

**Formulation:**

Interleukin -7 His is supplied in 1x PBS and 50% glycerol.

**Purity:**

Greater than 95.0% as determined by SDS-PAGE.

**Stability:**

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Please avoid freeze thaw cycles.